

Blood lipid and blood glucose test strip (dry chemical method) Instructions for use

Product Name

Common name: blood lipid and blood glucose test strip (dry chemical method)

Package Specification

Model: DiaCard-1 (TC + HDL-C + TG + GLU), Specification: 10 pcs/box, 20 pcs/box, 50 pcs/box

Model: DiaCard-2 (TC + HDL-C + TG), Specification: 10 pcs/box, 20 pcs/box, 50 pcs/box

Model: DiaCard-3 (TC + GLU), Specification: 10 pcs/box, 20 pcs/box, 50 pcs/box

Intended Use

It is used for in vitro quantitative detection of TC (total cholesterol), HDL-C (high density lipoprotein cholesterol), TG (triglycerides) and GLU (glucose) in human whole blood or serum samples. It is not used for glucose self-testing.

As one of the indexes of lipid metabolism, total cholesterol, high-density lipoprotein cholesterol and triglyceride are mainly used for the analysis of cardiovascular diseases. The test results can help patients find problems and seek medical treatment in time. Glucose can quickly detect the changes of glucose in the blood of diabetic patients, which has important guiding significance for life rules, activities, sports, diet and rational drug use, and helps patients find problems and seek medical treatment in time.

At present, the clinical and laboratory methods for detecting TC mainly include oxidase method, enzyme-coupled colorimetric method, COD-PAP method; the methods for detecting HDL-C mainly include HDL-CHO method, enzymatic method, and direct method; the methods for detecting TG mainly include enzymatic colorimetric method, end-point method, oxidase method; the methods for detecting GLU mainly include glucose oxidase method, dry chemical method.




Test Principle

After the blood sample is added to the test strip, the blood cells are filtered out and the serum diffuses to the reaction





layer. The substance to be measured in the serum reacts with enzymes and chemicals in the reaction layer and produces color changes. The color changes produced is proportional to the concentration. Then the analyzer uses the reflection method to read the results and calculates the content of the substance to be measured.

The reaction principle of each index is as follows:

Total cholesterol:





Triglycerides:



High-density lipoprotein cholesterol:

After precipitating LDL-C and VLDL-C with phosphotungstic acid in the presence of magnesium ions, HDL-C is measured using the same principle as the measurement of TC.

Glucose:



Composition

1. The blood lipid and blood glucose test strip contains the following active ingredients:

Name	Content	Name	Content
Cholesterol esterase	1.6U	Peroxidase	13U
Cholesterol oxidase	1.0U	Phosphotungstic acid	0.3mg
4-aminoantipyrine	150µg	Lipoprotein lipase	4.2U
N, N-disubstituted aniline	136µg	Glycerol kinase	2.4U
Glycerophosphate oxidase	1.7U	Adenosine triphosphate	50µg
Glucose oxidase	3.2U		

2. dial card

The dial card contains specific batch number of test strip for checking the test strip batch number.

Storage Conditions and Expiry Date

1. The test strip should be sealed and stored at 2°C ~ 30°C, and its period of validity is 12 months.
2. Test strip must be stored in the original packaging. Use the test strip immediately after taking out from package.
3. Refer to the packaging or label for the manufacture date and expiration date.

Applicable Analyzer

The test strip is applicable to Compass2000-1, Compass2000-2, Compass2000-3 dry biochemical analyzer produced by Jiangsu Konsung Bio-Medical Science And Technology Co., Ltd.

Specimen Requirements

1. Applicable blood samples include whole blood and serum.
2. If the sample contains anticoagulant, only whole blood anticoagulated by heparin or EDTA can be used.
3. Whole blood samples can be stored for 30 minutes at room temperature 25°C.
4. Serum samples can be stored at 2°C ~ 8°C for 7 days, and it can be stored at -20°C for 31 days.

Test Method

1. Prepare the necessary items for operation: blood sample, dry biochemical analyzer, test strip of corresponding type and other items required for operation.
2. Read this instructions for use of the dry biochemical analyzer carefully and be familiar with the operation of the analyzer.
3. Check if the dial card is consistent with the information on the test strip label.
4. Insert the blood lipid and blood glucose test strip into the correct position of the analyzer. The test strip DiaCard-1 requires 45µl blood specimen, DiaCard-2 requires 45µl blood specimen, and DiaCard-3 requires 20µl blood specimen.
5. The analyzer starts the test, and the result is displayed on the analyzer screen after reacting at 37°C for 3 minutes.
6. This product uses the national standard material GBW(E)090998 to calibrate total cholesterol, high-density lipoprotein cholesterol and triglycerides; the national

standard material GBW(E)091004 to calibrate glucose.

7. This product uses Randox quality control serum (item number: HE1532 and HN1530) for quality control.

TG	0.51~7.34	r≥0.975
HDL-C	0.39~2.59	r≥0.975
GLU	2.0~18.0	r≥0.975

Reference Intervals

This product tests and analyzes the serum samples of healthy adults. The upper limit of the reference interval for TC is the 95th percentile; the lower limit of the reference interval for HDL-C is the 5th percentile; the upper limit of the reference interval for TG is the 95th percentile; the reference interval for glucose is adopted $x \pm 1.96s$, and the reference range of each item is as follows:

Test Items	Sample cases	Reference range
TC	155	<5.2mmol/L
HDL-C	155	>1.0mmol/L
TG	155	<1.7mmol/L
GLU	155	3.9mmol/L~6.1mmol/L

Each laboratory shall determine the applicability of the reference interval through tests, and establish its own reference interval range based on the patient population being tested when necessary.

Test Result

Interference factors:

- Bilirubin (>20mg/dL), vitamin C (>10mg/dL), uric acid (>10mg/dL) will make the test results low.
- Dopamine (>1.8mg/dL) will make the test results low.

Limitation of Test Method

Diagnosis and treatment should not rely solely on the test result. Clinical history and other laboratory tests should be considered.

Performance Specification

- Blank limit
All items shall not be higher than 0.3mmol/L.
- Linearity

Test item	Linearity Range (mmol/L)	Correlation requirements
TC	2.59~12.93	r≥0.975

3. Repeatability

- The CV (Coefficient of Variation) of total cholesterol is ≤7.5%;
- The CV (Coefficient of Variation) of triglycerides is ≤6%;
- The CV (Coefficient of Variation) of HDL-C is ≤7.5%;
- The CV (Coefficient of Variation) of glucose is ≤5%.

4. Inter-batch difference

- The inter-batch difference of total cholesterol is ≤15%;
- The inter-batch difference of triglyceride is ≤15%;
- The inter-batch difference of HDL-C is ≤15%;
- The inter-batch difference of glucose is ≤10%.

5. Accuracy

- The relative deviation of total cholesterol is ≤±10%;
- The relative deviation of triglycerides is ≤±15%;
- The relative deviation of HDL-C is ≤±15%;
- Relative deviation of glucose is ≤±10%.

Notes

- This product is for In Vitro Diagnostic use only.**
- Use the test strip prior to the expiration date.**
- During the operation, fill the test strip in one continuous process. Do not refill repeatedly. If the test strip cannot be filled with enough blood specimen in one continuous process, use another new test strip to perform the test again.**
- This product is for single-use only. Always handle blood specimens with care as they may be infectious. Consult local environment authorities for proper disposal. Always wear protective gloves when handling blood specimens and test strips with specimen.**

Symbol and Explanation

Symbol	Explanation	Symbol	Explanation
	In vitro diagnostic medical device		Date of manufacture
	Consult instructions for use		Manufacturer

	Caution		Batch code
	Temperature limit		Do not reuse
	Use-by date		UKCA mark
	Name and Address of UK responsible person		Authorized representative in the European community
	CE mark		

Bibliography

- Joint committee for guideline revision of Chinese Adult Dyslipidemia Prevention, Chinese guidelines for the management of dyslipidemia in adults (2016 Revised Edition). Chinese Circulation Journal, 2016, 31: 937-953.
- Friedewald et al. Clin Chem. 1972, 8(6):499-502.



Jiangsu Konsung Bio-Medical Science And Technology Co., Ltd.

NO.8, Shengchang West Road, Danyang Development Zone, Jiangsu Province, 212300, P.R.China



SUNGO Certification Company Limited

3rd floor, 70 Gracechurch Street, London. EC3V 0HR



Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80,20537 Hamburg,Germany



Fitech UK Ltd

About this Instructions for use

P/N: 138040411

Version: 1.0

First edition/Release date: June of 2024